




ORYGEN DUAL VALVE

Πρωτοποριακή συσκευή αναπνευστικής
φυσιοθεραπείας

- Η συσκευή ORYGEN DUAL VALVE είναι κατάλληλη για την εκπαίδευση των αναπνευστικών μυών με μία αναπνοή
- Είναι συσκευή 2 σε 1
- Διαθέτει το πλεονέκτημα εκπαίδευσης των εισπνευστικών μυών καθώς και των εκπνευστικών
- Δυνατότητα ρύθμισης αντιστάσεων
Για εισπνευστικές έως 90 cmH₂O και εκπνευστικές έως 60 cmH₂O
- Καταγραφή χρόνου εξάσκησης με ημερήσιο πρόγραμμα (training diary)
- Χορηγείται κατόπιν ιατρικής εντολής



ORYGEN DUAL VALVE

LUNG TRAINER

SEMANA WEEK	DOMINGO / SUNDAY				LUNES / MONDAY				MARTES / TUESDAY			
	TIME	MIN	TIME	MIN	TIME	MIN	TIME	MIN	TIME	MIN	TIME	MIN
1												
2												
3												
4												
5												
6												
7												
8												
9												
10												
11												
12												
13												
14												
15												
16												
17												
18												
19												
20												



GENERAL DESCRIPTION

THE ORYGEN DUAL VALVE SYSTEM IS AN OPEN CIRCUIT WITH INSPIRATORY AND EXPIRATORY PORTS, FOR INSPIRATORY AND EXPIRATORY MUSCLES TRAINING.

ORYGEN Dual Valve allows for imposing external pressure loads on both the inspiratory and expiratory circuit. The system has been designed to act independently in both inspiratory and expiratory phases. The imposed loads can vary according to the functional characteristics of the muscle groups in each patient. The valve is factory adjusted at a minimal threshold load.

SPECIFICATION

Maximum values of the pressure imposed by the system:

- inspiratory circuit: 90 cmH₂O
- expiratory circuit: 60 cmH₂O

Regular intervals of both cases: 10 cm/H₂O