



*FitMax*<sup>™</sup> CPAP total face mask is the best alternative for those patients who are not able to obtain a good comfort with nasal or full-face masks, or those with sensitive skins on face or nasal area or who is claustrophobia.

The appropriate perimeter of face reduces the pressure sensitivity and enhances facial contour with comfort.

# *FitMax*<sup>™</sup>

- Medical grade silicone sealing cushion provides optimal comfort for the patients.
- Superlight weight and clear view to ease anxiety and claustrophobia.
- Slide-in clips for headstrap design for easy-to-use.
- **Anti-Asphyxia Elbow** is for patient safety; optional **Standard Elbow** good for ventilator application.
- White and blue double faced soft sponge bonnet fits infants' demand.
- Available sizes: XS, S, M, L, XL



## Feature

Comfortable Headstrap

Washable design

Medical Grade Silicone Cushion

Swivel connector

Large Clear Surface

Pressure Measurement Port

Patented DE, TW, CN

Slide-in Clip for easy use.

DE Pat. 202014104150U1


TW Pat. M490327

ZN Pat. ZL201420181488.8

















## Circuits










With a choice of the **Standard Elbow** for dual-limb systems and the **Anti-Asphyxia Elbow** for single-limb systems, clinicians can use the **FitMax™** Mask with any non-invasive ventilator.



## FitMax™ CPAP Total Face Mask - Single-Patient Use

| Component |    | FitMax™<br>CPAP MASK   | Connector   |   |   | Head Strap |
|-----------|----|--|---|---|---|------------|
|           |    |  | CA-63612<br>Standard Elbow  | 22M/15F, 15M<br>15F/22M, 15F/22M  | CA-63611<br>Anti-Asphyxia Elbow   |            |
| Item      |    |             |  |  |  |            |
| Infant    | XS |  CA-63603   | •   | •   | •   | •          |
|           |    |  CA-63605   | •   |   | •   | •          |
| Child     | S  |  CA-63607   | •   | •   | •   | •          |
|           |    |  CA-63609   | •   |   | •   | •          |
| Adult     | M  |  CA-63600   | •   | •   |   | •          |
|           |    |  CA-63601   | •   |   | •   | •          |
|           | L  |  CA-63602  | •   | •   |   | •          |
|           |    |  CA-63604 | •   |   | •   | •          |
|           |    |  CA-63606 | •   | •   |   | •          |
|           | XL |  CA-63608 | •   |   | •   | •          |

## FitMax™ Reusable CPAP Total Face Mask

| Component |    | FitMax™<br>CPAP MASK   | Connector   |   |   | Head Strap |
|-----------|----|--|---|---|---|------------|
|           |    |  | CA-63612<br>Standard Elbow  | 22M/15F, 15M<br>15F/22M, 15F/22M  | CA-63611<br>Anti-Asphyxia Elbow   |            |
| Item      |    |           |  |  |  |            |
| Infant    | XS |  CA-63621 | •   | •   | •   | •          |
| Child     | S  |  CA-63622 | •   | •   | •   | •          |
| Adult     | M  |  CA-63623 | •   | •   | •   | •          |
|           | L  |  CA-63624 | •   | •   | •   | •          |
|           | XL |  CA-63625 | •   | •   | •   | •          |

## Accessories

### Headstrap



**HE-74110** For Infant (XS)



**HE-74210** For Pediatric (S)



**HE-74310** For Adult (M, L, XL)

### Connector



**CA-63611** Anti-Asphyxia Elbow



**CA-63612** Standard Elbow


### Applicable Connector



**CT-14316** Swivel Straight Connector



**CT-14317** Straight Connector



### Durable CPAP Tubing W/TPR Cuffs



ID: 22mm

**CT-10825** 1.2m

**CT-10827** 1.5m

**CT-10900** 1.8m

**CT-10830** 2.0m



#### **BESMED HEALTH BUSINESS CORP.**

No. 2, Lane 106, Wu-Kong 3rd Rd, Wu-Ku District, New Taipei City, Taiwan. 24888

**Tel:** +886-2-2290-3959

**Fax:** +886-2-2299-9076

**E-mail:** [info@besmed.com](mailto:info@besmed.com)

**Web-site:** [www.besmed.com](http://www.besmed.com)



2018-07 Rev.3