

Your cardio-respiratory partner

Data you can rely on. People you can trust.



# **Vitalograph Peak Flow Meters**

## Peak flow monitoring for adults and children

Vitalograph peak flow meters combine high quality and attractive styling with great value. These robust and easy to clean peak flow meters conform to the ISO standard for Peak Expiratory Flow Meters.

## myPEF – the only web enabled peak flow meter

myPEF™ is supplied with a unique ID card allowing users free access to a secure website where they can enter their PEF readings, eliminating the need for paper records and supporting data accuracy and effective asthma management. Suitable for use with any web enabled device.

- Easy to read graphical data display featuring data trending, personal best and asthma management plan
- Report generation featuring plotting, trending and time stamping
- Olinicians can be granted access to view patient's PEF results online
- Precise 100mm scale longer than any other PFM's scale
- Integral dust cover for hygiene and protection
- Lightweight, pocket size design for convenient portability

#### asmaPLAN

asmaPLAN™ is the first and original home use personal monitor for asthma management. Developed in the early 1970's it has evolved to become the gold standard of simple mechanical peak flow meters. The lightweight, proven and truly robust asmaPLAN is highly accurate and easy to use by both adults and children.

- Proven
- Low cost
- Suitable for adults and children
- Easy to use
- Highly accurate
- Robust







### Technical specification:

**Product:** Peak Flow Meter **Model Number:** 4300

**Parameters Measured:** PEF (L/min)

**Device Dimensions:** 

138mm x 56mm x 23mm (myPEF) 152mm x 61mm x 36mm (asmaPLAN)

Nett Weight: 79g (myPEF); 96g (asmaPLAN)

**Environmental Data:** Storage humidity: 10% to 95%, Storage temperature: 0-50°C

**Recommended Operating Temperature Range:** 17 - 37°C; Design limits 10-40°C

**Ambient Humidity:** 0% to 99%

**Medical Safety Standard:** Medical Devices

Directive 2007/47/EC

**Designed and Manufactured Under:** ISO 13485:2003, FDA 21 CFR 820, CMDR SOR 98/282 and JPAL

#### Ordering info:

43950 4300 Peak Flow Meter myPEF Web Enabled

43602 4300 Peak Flow Meter asmaPLAN

43703 4300 Peak Flow Meter asmaPLAN Child

**43982** 4300 Peak Flow Meter myPEF Disposable (100)

**20242** 2024 SafeTway® Mouthpieces (200)

**29805** Mouthpiece Adaptors myPEF (96)

Vitalograph is a world leading provider of outstanding quality cardio-respiratory diagnostic devices, clinical trial services and medical equipment servicing.

With a pioneering heritage of excellence spanning half a century Vitalograph continues to make valuable contributions to effective medical care and enhanced quality of life.



Your cardio-respiratory partner

Vitalograph Ltd. UK Maids Moreton, Buckingham

MK18 1SW England

Phone: (01280) 827110 Fax: (01280) 823302

e-mail: sales@vitalograph.co.uk www.vitalograph.co.uk

Vitalograph Inc. 13310 West 99th Street Lenexa, Kansas, 66215

Toll Free: 800 255 6626 Phone: (913) 730 3200 Fax: (913) 730 3232

e-mail: vitcs@vitalograph.com www.vitalograph.com Vitalograph Ltd. Export Maids Moreton, Buckingham

MK18 1SW England

Phone: +44 1280 827120 Fax: +44 1280 823302 e-mail: sales@vitalograph.co.uk www.vitalograph.eu

Vitalograph (Ireland) Ltd. Gort Road Business Park Ennis, Co. Clare

Ireland Phone: +353 65 6864100 Fax: +353 65 6829289 e-mail: sales@vitalograph.ie www.vitalograph.ie Vitalograph GmbH Rellinger Straße 64a D-20257 Hamburg Germany Phone: (040) 547391-0 Fax: (040) 547391-40 e-mail: info@vitalograph.de www.vitalograph.de

Vitalograph Ltd. Hong Kong Unit 230, Corporation Park 11 On Lai Street, Shatin Hong Kong Phone: +852 2117 2678 Fax: +852 2117 2679 e-mail: sales@vitalograph.hk

www.vitalograph.hk